Examples showing different grips and positions available in advanced handprostheses and in what simple tasks they may be used.

Grip/position		Task/objects
Power grip		Carry a glass, handle a pepper mill.
Tripod pinch	0	Pull up zippers, play cards.
Pinch grip precision	8	Fine manipulation of objects. Hold a needle.
Precision grip hand open	10	Fine manipulation of objects. Pick up coins from a purse.
Precision grip hand closed	0	Fine manipulation of objects. Pick up small objects from a desk or a counter, e.g. coin, paperclip.
Key pinch		Eat with spoon, carry paper, hold a credit card or a key.
Pointing indexfinger		Type on keyboard, press a button, e.g. credit card code.
Active index		Handle spray bottle, power tool or hair dryer.
Computer mouse grip		Use a computer mouse.
Column position		Dress the upper limb – pull into a sleeve. Push larger buttons, e.g. door openers.
Hook grip	D	Carry a bag or briefcase.
Adduction grip		Hold small objects, cards or cutlery.
Open palm position	Y	Carry plates, cell phone or wallet. Stabilize e.g. bread loaf when cutting into slices.
Relaxed hand position	è	Walk in a natural and relaxed way.
Gestures	b	Express your feelings in different ways.